

Ice rack Ice On




Ice-On is a tube stand which cools samples quickly compared to placing onto ice directly. It reduces contamination risk compared with tube samples placed directly on ice.

More information
https://www.funakoshi.co.jp/exports_contents/80502



For FACS Mini Cell Strainers



Mini Cell Strainers are small size cell strainers, which fit for 1.5 ml/2.0 ml tubes, 5 ml sampling tubes and 15 mL conical tubes. It facilitates speedy cell separation from primary tissues or cell clumps.

More information
https://www.funakoshi.co.jp/exports_contents/80727



Grab for something hot Hot Hand



Your finger is protected when grabbing something hot. Just insert your fingers into the pockets of Hot Hand and bend like left picture.

6 sizes are available.

More information
https://www.funakoshi.co.jp/exports_contents/46146



NOTE

※ All products here are research use only, not for diagnostic use.
 ※ Specs might be changed for improvement without notice.
 ※ Numbers after “#” represents product code.

※ Company name and product name are trademark or registered mark.
 ※ Please contact your local distributors for orders, quote request and inquiry.

Your Local Distributor

Bio-REV Pte. Ltd.

Block 211 Woodlands Avenue 9
 #08-78/79 Woodlands Spectrum II
 Singapore 738 960
 TEL : +65-6555-9001
 FAX : +65-6555-9030
 Email: sales@bio-rev.com
 Technical Support: techserv@bio-rev.com

Funakoshi Co., Ltd. (JAPAN)

Phone: +81-3-5684-6296
 Mail: export@funakoshi.co.jp



\Contact/



\Newsletter Signup/

P4

FUN-2009-F06 (2020.09)

funakoshi
 FRONTIERS IN LIFE SCIENCE

www.funakoshi.co.jp/exports
 Research Use Only



For your comfortable lab routine

Your Lab Essentials

How is your lab life? Do you feel comfortable with your tools which you routinely use? We, Funakoshi, help your routine task to be easy, smooth and comfortable. Let us introduce when our items support your research with “My Lab Routine”!

My Lab Routine

Funakoshi items

Morning 	Lab clean	- Autoclave bags that decrease the smell P2
	See my cells	- Plastic Pasteur pipettes P2 - Stretchy covering film P2
	Start WB	- Prepare reagents without frustration P2 - Easy to use! Protein marker of SDS-PAGE P3
Afternoon 	Lunch break	
	Start IHC	- Serial section floating bath P3 - Let's share a photomicrograph with mentor P3
	Take a rest	
Evening 	Start PCR	- 10 minutes! Rapid and Easy DNA Extraction P3 - Keep samples cool P4
	Read paper	
	Lab finish	

Other items

- For FACS, Mini Cell Strainers **P4**
- Grasp something hot safely **P4**

P1

Lab clean

Autoclave bags with special effect



There are two special types.

One is **autoclave bag with deodorant component in the film** which decreases the smell of waste.


Another is **autoclave double bag** which is solidly bag with PP(Polypropylene) and PE(Polyethylene).

More information

https://www.funakoshi.co.jp/exports_contents/81429



Protein marker



Protein MultiColor Stable II

Pre-stained protein ladder for broad range of proteins (8 - 230 kDa), stable in refrigerator (at 4°C) for 1 year !

More information

https://www.funakoshi.co.jp/exports_contents/81388



Cell culture

Plastic Pasteur pipettes



Plastic Pasteur pipettes are suitable for cell culture, vacuuming liquids, transfer small liquid quantities and other applications.

More information

https://www.funakoshi.co.jp/exports_contents/81411



Section item

Serial section floating bath



Serial section floating bath is useful bath (vat) for floating serial sections in order of sectioning. This has appropriate depth for easy collecting of sections.

More information

https://www.funakoshi.co.jp/exports_contents/81510



Stretchy covering film

Labopita



Labopita is a new autohesion sealing film, which has a great capability to chemical resistance, moisture-proof property, sealing property and stretch property.

More information

https://www.funakoshi.co.jp/exports_contents/81178




Microscope imaging by iPhone

i-NTER LENS



i-NTER LENS (interlens) is an innovative attachment for taking microscope images by iPhone. It enables to share the photomicrograph.



More information

https://www.funakoshi.co.jp/exports_contents/81175



Spatel / Tube / Weighing dish

Antistatic Consumables



(A) Antistatic spatel
(B) General spatel (polypropylene)
(C) General spatel (stainless)

Our **spatel, tube and weighing dish** contains high molecular type antistatic additive, then ideal for dispensing, weighing and storing of microparticles. Compared to low molecular type antistatic additive,

antistatic effect is kept semipermanently.


More information

https://www.funakoshi.co.jp/exports_contents/80984



DNA extraction

Easy DNA Extraction Kit (Version 2)



P3 DNA extraction from mouse tail sample

Kaneka Easy DNA Extraction Kit (Version 2) enables you to extract template DNA which can be used for downstream applications such as PCR or Real time PCR. Just two components and simple protocol makes your DNA extraction process more **rapid (around 10 minutes)**.

More information

https://www.funakoshi.co.jp/exports_contents/80099

