

Steritainer Plus

Gamma Beam Sterilized

For more information : https://www.funakoshi.co.jp/exports_contents/81430

3D Single Use Sterile Bag is suitable for the storage / transportation of liquids which requires sterile environment in a container.

Features

-Space saving

Compact package individually wrapped.
Stores flat saving space when unopened.

-High-Strength Design

Self-supporting design, stirring can be carried out easily without transferring the solution to glass containers.

- Easy to Use

Special embossed coating prevents fold wrinkles and scratches.
A transparent window allows visualization of solution.



<Option>

3-Port Connector / 3-Port Connector ASSY

-For discharge solution and air in closed sterile system.

Class	Port OD and No. of Ports
1/4 inch Type	1/4 inch Type x 3
3/8 inch Type	3/8 inch Type x 2
	1/4 inch Type x 1

3-Port Connector



3-Port Connector ASSY



Product Information

[Manufacturer : SEC]

Product Name	Code	Capacity(ℓ)	Container Inner	Size	Packing Dimension (mm)
			Dimension(mm)		
Steritainer Plus	SR-05S	5	190 × 190 × 190H	10Sheets	Outer Box:400 × 295 × 295H Inner Box:279 × 279 × 183H
	SR-10S	10	230 × 230 × 230H		
	SR-20S	20	285 × 285 × 285H		
Steritainer Plus 3ports connector included	SR-05P	5	190 × 190 × 190H		
	SR-10P	10	230 × 230 × 230H		
	SR-20P	20	285 × 285 × 285H		
Steritainer Plus 3ports connector ASSY (1/4inch) included	SR-05A(1/4)	5	190 × 190 × 190H	5Sheets + 5ASSY	
	SR-10A(1/4)	10	230 × 230 × 230H		
	SR-20A(1/4)	20	285 × 285 × 285H		
Steritainer Plus 3ports connector ASSY (3/8inch) included	SR-05A(3/8)	5	190 × 190 × 190H		
	SR-10A(3/8)	10	230 × 230 × 230H		
	SR-20A(3/8)	20	285 × 285 × 285H		

Autoclave bags with special effect

For more information : https://www.funakoshi.co.jp/exports_contents/81429

Have a comfortable laboratory life!

Easy packing !

Closure type: Zip tie



**There are two
special types.**

Type1. Autoclave bag with deodorant component in the film

Decrease the smell of waste.

S-Clinica Pack M

Size: 50 x 70 cm

Thickness: 0.08mm

Material: PP(Polypropylene)

Autoclave Temperature : up to 130 °C

S-Clinica Pack L

Size: 62 x 82 cm

Thickness: 0.08mm

Material: PP(Polypropylene)

Autoclave Temperature : up to 130 °C

Type2. Autoclave double bag

Double layer solidly bag with PP(Polypropylene)
and PE(Polyethylene).

*No deodorant component contained.

S-Clinica Pack WM

Size: Outer 50 x 70 cm , Inner 49.5 x 72 cm

Thickness: Outer 0.06mm , Inner 0.04mm

Material: Outer PP(Polypropylene) , Inner PE(Polyethylene).

Autoclave Temperature : up to 130 °C

S-Clinica Pack WL

Size: Outer 62 x 82 cm , Inner 61.5 x 84 cm

Thickness: Outer 0.06mm , Inner 0.04mm

Material: Outer PP(Polypropylene) , Inner PE(Polyethylene).

Autoclave Temperature : up to 130 °C

Product Information

[Manufacturer : SEC]

Product Name	Code	Size	Storage
S-Clinica Pack M	C-CP11	200 /case	RT
S-Clinica Pack L	C-CP12	200 /case	
S-Clinica Pack WM	C-CE21	100 /case	
S-Clinica Pack WL	C-CE22	100 /case	

NOTES

- ※ All products here are research use only, not for diagnostic use.
- ※ Specs might be changed for improvement without notice.
- ※ Numbers after “#” represents product code.
- ※ Company name and product name are trademark or registered mark.
- ※ Please contact your local distributors for orders, quote request and inquiry.

Your Local Distributor



Bio-REV Pte. Ltd.

Block 211 Woodlands Avenue 9

#08-78/79 Woodlands Spectrum II

Singapore 738 960

TEL : +65-6555-9001

FAX : +65-6555-9030

Email: sales@bio-rev.com

Technical Support: techserv@bio-rev.com

Funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,

Tokyo 113-0033 JAPAN

Phone : +81-3-5684-6296

Fax : +81-3-5684-6297

Email : export@funakoshi.co.jp

SEC-2002-F08 (2020.2)